

# CANANOLAB 1.5.1

## *Installation Guide*



May 20, 2010

# Table of Contents

<b>INTRODUCTION TO CANANOLAB .....</b>	<b>1</b>
Targeted Developer.....	1
General System Requirements .....	2
Obtaining the caNanoLab Source .....	2
 <b>DATABASE TECHNOLOGY .....</b>	 <b>3</b>
Initializing the MySQL Database .....	3
Follow these steps to initialize your MySQL database system:.....	3
Data Migration .....	4
 <b>CANANOLAB WEB APPLICATION.....</b>	 <b>5</b>
User Provisioning Tool (UPT).....	9
 <b>CANANOLAB GRID DATA SERVICE.....</b>	 <b>13</b>
 <b>CONTACTING APPLICATION SUPPORT .....</b>	 <b>16</b>

## Introduction to caNanoLab

Welcome to the cancer Nanotechnology Laboratory (caNanoLab) 1.5.1 Installation Guide. caNanoLab is a data sharing portal designed to facilitate information sharing in the biomedical nanotechnology research community to expedite and validate the use of nanotechnology in biomedicine. caNanoLab allows researchers to share information on nanomaterials including the composition of the nanomaterial, the functions (e.g. therapeutic, targeting, diagnostic imaging) of the particle, the characterizations of the nanomaterial from physico-chemical (e.g. size, molecular weight, surface) and in vitro (e.g. cytotoxicity, blood contact) assays, and the protocols of these characterization.

As of release 1.1, the caNanoLab domain model has been caGrid enabled. In other words, a caNanoLab grid data service can be deployed and registered with caGrid production index server, allowing sharing of public nanomaterial information across the caGrid. In release 1.5, the caNanoLab grid service has been updated to use caGrid 1.3. For more information, see <https://cabig.nci.nih.gov/workspaces/Architecture/caGrid>.

---

## Targeted Developer

The caNanoLab application development is best suited for an experienced Java developer who has some familiarity with the following J2EE and related technologies:



### Important Background Knowledge

- **Unix/Linux environment, Windows XP environment or Mac OS environment (Configuring environment variables; Installing Ant, JDK, Apache Tomcat and JBoss servers)**
  - **Ant build scripts**
  - **J2EE web application development using the Struts framework, Servlet/JSP's, JavaScript**
  - **J2EE middle-ware technologies such as n-tier service-oriented architecture and software design patterns**
  - **Hibernate Java persistence framework**
  - **caGrid 1.3 infrastructure (for understanding how a caNanoLab grid data service functions)**
-

## General System Requirements

The following open source technologies power a caNanoLab web application:

### Open Source Technologies

- Java Software Development Kit (JDK) version 5.0  
[http://java.sun.com/javase/downloads/index\\_jdk5.jsp](http://java.sun.com/javase/downloads/index_jdk5.jsp)
  - JBoss version 4.0.5 <http://labs.jboss.com/jbossas/downloads>
  - Jakarta Ant version 1.6.5 <http://archive.apache.org/dist/ant/binaries/>
  - MySQL version 5.0.x <http://dev.mysql.com/downloads/mysql/5.0.html>
- 

The caNanoLab web application has been tested within NCICBIIT against JBoss servers (version 4.0.5) hosted on Windows XP and RedHat Linux systems, and against MySQL 5.0.45 databases hosted on RedHat Linux systems and MySQL 5.0.45 on Windows XP systems. Prior to release 1.2.1, the caNanoLab web application had been tested against Oracle 9i databases hosted on Sun Solaris systems, and Oracle 10g XE database hosted on Windows XP.

Download each of the tools listed in the bulleted list above and follow the installation instructions provided with each respective product for your environment. Assistance from a MySQL database administrator is expected to properly configure the MySQL database. For MySQL database configuration and maintenance, it is also helpful to download the MySQL 5.0 suite of GUI tools at <http://dev.mysql.com/downloads/gui-tools/5.0.html>.

---

### Grid Service

Setting up a caNanoLab grid service is optional. It is recommended that you install the grid service only after you successfully install the web application and become familiar with submitting and searching data through the web application.

---

### Grid Service Downloads

In release 1.5.1, there are no separate technology downloads required for installing the caNanoLab grid service, and most have been packaged within the caNanoLab source distribution file.

---

## Obtaining the caNanoLab Source

### Source Code

The caNanoLab web application source code and the caNanoLab grid service deployment war file are distributed as a zip file named `caNanoLab_1.5.1.zip` from the NCICBIIT GForge site at [http://gforge.nci.nih.gov/frs/?group\\_id=69](http://gforge.nci.nih.gov/frs/?group_id=69) under release 1.5

---

## Database Technology

### Assumptions and Requirements

The caNanoLab source distribution `caNanoLab_1.5.1.zip` has been downloaded from the NCICBIIT GForge site. A MySQL database has been set up on a system (local or remote) with an administrator (or root) account. You should be familiar with the process of creating InnoDB tablespaces, databases creating users and assigning privileges to users.

### For a Previous Installation

If you have installed caNanoLab prior to release 1.5.1 against an MySQL database, or release 1.2 or release 1.1.1 or release 1.1 against an Oracle database, and have associated production data in these schemas and you would like to continue to use the same data for caNanoLab release 1.5.1 in MySQL, review the following database initialization steps, then go directly to *Data Migration* on page 4. If you have a previously installed caNanoLab application running on a JBoss server, stop the server before running the database scripts.

### For a New Installation

If you are installing the caNanoLab application for the first time or want to install a new schema for release 1.5.1, follow the steps below to set up the required MySQL schema objects and the seed data for release 1.5.1.

## Initializing the MySQL Database

Follow these steps to initialize your MySQL database system:

Step	Action
1	<p>Extract the <code>caNanoLab_1.5.1.zip</code> to a location on your local system, for example, <code>C:\caNanoLab_1.5.1</code>. This location is referred as <code>&lt;CANANOLAB_SOURCE&gt;</code> throughout the document. Verify that the following seven SQL scripts exist in the directory <code>&lt;CANANOLAB_SOURCE&gt;/db/1.5.1</code>:</p> <ul style="list-style-type: none"> <li><code>mysql_seed.sql</code></li> <li><code>create_database_and_user.sql</code></li> <li><code>app_schema.sql</code></li> <li><code>insert_common_lookup.sql</code></li> <li><code>csm_schema.sql</code></li> <li><code>csm_priming_data.sql</code></li> <li><code>app_csm_priming_data.sql</code></li> </ul> <p><b>Note:</b> The scripts <code>create_database_and_user.sql</code> create a new 1.5.1 schema <b>canano</b> and a new database user <b>cananolab_app</b> with required privileges and a password <b>go!234</b>. If you'd like to create a different database user, be sure to modify this script and make a note of the information for later use</p>

Step	Action
	<p>in setting up the caNanoLab application.</p> <p>User administration and schema privileges can also be managed by the MySQL Administrator GUI as a part of the MySQL 5.0 suite of GUI tools.</p>
2	<p>Log into the MySQL database as the database administrator (root) and go to the directory <code>&lt;CANANOLAB_SOURCE&gt;/db/1.5.1</code> and execute the script <code>mysql_seed.sql</code> to set up a schema named <b>canano</b> and its required seed data.</p> <p><i>Example:</i> On a Windows XP system hosting a MySQL 5.0.45 database, you would issue the following commands at the DOS prompt (assuming the <b>root</b> account has password <b>rootpass</b>):</p> <pre>C:\&gt;cd C:\caNanoLab_1.5.1\db\1.5.1</pre> <pre>C:\caNanoLab_1.5.1\db\1.5.1&gt;mysql -h localhost -u root -prootpass &lt;mysql_seed.sql</pre>

## Verification

Once the MySQL database has been created, either through new setup or through data migration (described in the next section), verify that the following numbers of database objects are created:

- Tables 67

*Example:* Issue the following query at the MySQL prompt, logging in as **cananolab\_app**:

```
mysql> select count(*) from information_schema.tables
where table_schema='canano' and table_type='BASE TABLE';
```

## Data Migration

### NOTE:



If you are installing caNanoLab 1.5.1 for the first time or installing a new caNanoLab release 1.5.1 schema, you can skip this section.

This installation guide only discusses the steps for migrating from release 1.5 in MySQL to release 1.5.1 in MySQL. If you have previously installed caNanoLab in MySQL prior to release 1.5, you'd have to migrate to release 1.5 first, one release at a time. Please see the install instructions for each release at [http://gforge.nci.nih.gov/frs/?group\\_id=69](http://gforge.nci.nih.gov/frs/?group_id=69) for details.

Follow these steps to complete the required data migration from release 1.5 in MySQL to release 1.5.1 in MySQL:

Step	Action
1	<p>Extract the caNanoLab_1.5.1.zip to a location on your local system, for example, C:\caNanoLab_1.5.1. This location is referred as &lt;CANANOLAB_SOURCE&gt; throughout the document. Verify that the following SQL scripts exists in the directory &lt;CANANOLAB_SOURCE&gt;/db/1.5.1/data_migration:</p> <ul style="list-style-type: none"> <li>from_1.5.sql</li> </ul>
2	<p>Log into the MySQL database as the database administrator (root) and go to the directory &lt;CANANOLAB_SOURCE&gt;/db/1.5.1/data_migration and execute the following script in the order specified:</p> <ul style="list-style-type: none"> <li>from_1.5.sql</li> </ul> <p><i>Example:</i> On a Windows XP system hosting a MySQL 5.0.45 database, you would issue the following commands at the DOS prompt (assuming the <b>root</b> account has password <b>rootpass</b>):</p> <pre>C:\&gt;cd C:\caNanoLab_1.5.1\db\1.5.1\data_migration C:\caNanoLab_1.5.1\db\1.5.1\data_migration&gt;mysql -h localhost -u root -prootpass &lt;from_1.5.sql</pre>

After data migration, refer to the Verification section on page 4 to verify that the migration has been successful. There will be 68 instead of 67 schema objects if had migrated 1.5 from 1.4.1.2 before. The extra table `instrument_to_review` contains instruments that were not automatically migrated and requires manual review in release 1.5.

## caNanoLab Web Application

**Assumptions  
and  
Requirements**

The database has been installed and verified, as described on pages 3 through 5. Ant has been installed. The JBoss application server has been installed on a system (local or remote) and can be started at a designated port. JBoss install directory is referred as `<JBOSS_HOME>` in the document. The JBoss application server host URL is referred as `<APP_SERVER_URL>`, for example, `cananolab-dev.nci.nih.gov:19080`. It is assumed that the default configuration is used for deploying caNanoLab web archive files. For example, in JBoss 4.0.5, the default configuration is located at the directory `<JBOSS_HOME>/server/default`. The caNanoLab web archive file shall be deployed at the directory `<JBOSS_HOME>server/default/deploy`.

It is possible to configure Apache server to interface with the JBoss server and set up a virtual host for the caNanoLab application, if you need assistance, please contact NCICBIIT Application Support (info provided at the end).

**Installation and  
Deployment**

Follow these steps to install and deploy caNanoLab.

Step	Action
1	<p>Set up an environment variable <code>JBOSS_HOME</code> to point to the JBoss installation directory.</p> <p><b>Note:</b> Either double back slashes <code>\\</code> or a single forward slash <code>/</code> should be used as the file separator if working on Windows platform, for example, <code>JBOSS_HOME=C:/jboss-4.0.5.GA</code> (<code>C:\jboss-4.0.5.GA</code> would fail).</p>
2	<p>Edit the Ant build properties file <code>build.properties</code> at <code>&lt;CANANOLAB_SOURCE&gt;</code> by specifying values for the following three properties:</p> <ul style="list-style-type: none"> <li>a. <code>file.repository.dir</code>: a directory on the system that hosts the JBoss application server for storing uploaded files, for example, <code>C:/caNanoLab/fileupload</code>.  <p><b>Note:</b> This directory should be writable by the user that starts the JBoss server, and this directory should be created prior to starting the application. Either double back slashes <code>\\</code> or a single forward slash <code>/</code> should be used as the file separator if working on Windows platform.</p> </li> <li>b. <code>application.owner</code>: the owner of the caNanoLab release 1.5.1 installation instance, for example, <code>NCICBIIT</code>.</li> <li>c. <code>admin.email</code>: the email address that receives requests for new user accounts, for example, <code>NCICBIIT application support email for the NCICBIIT installation</code>.</li> </ul>



Step	Action
3	<p>Edit the Ant build properties file <code>deploy.properties</code> at <code>&lt;CANANOLAB_SOURCE&gt;</code> by specifying values for the following four properties:</p> <ol style="list-style-type: none"> <li><code>database.host</code>: the host name of the server hosting the database, for example, <code>localhost</code></li> <li><code>database.port</code>: the database port number, for example, <code>3306</code></li> <li><code>database.user</code>: the database user name used in the caNanoLab application, for example, <code>cananolab_app</code>. Please refer to page 3.</li> <li><code>database.password</code>: the password for the database user specified above, for example, <code>go!234</code></li> </ol>
4	<p>Execute the Ant build script <code>build.xml</code> located at <code>&lt;CANANOLAB_SOURCE&gt;</code> with the target <code>deploy-application</code>.</p> <p><i>Example:</i> Issue the following commands to execute the Ant script:</p> <pre>C:\&gt;cd C:\caNanoLab_1.5.1 C:\caNanoLab_1.5&gt;ant deploy-application</pre> <p>Successful execution of the Ant script creates two deployable web archive war files in the directory <code>&lt;CANANOLAB_SOURCE&gt;/target</code>: <code>caNanoLab.war</code> and <code>upt.war</code>, and places these two wars to the JBoss deploy directory, for example, <code>&lt;JBOSS_HOME&gt;/server/default/deploy</code>.</p>
5	<p>Execute the Ant build script <code>build.xml</code> located at <code>&lt;CANANOLAB_SOURCE&gt;</code> with the target <code>prepare-csm-conf</code>.</p> <p><i>Example:</i> Issue the following commands to execute the Ant script:</p> <pre>C:\&gt;cd C:\caNanoLab_1.5.1 C:\caNanoLab_1.5&gt;ant prepare-csm-conf</pre> <p>Successful execution of the Ant script sets up the required CSM 4.1 configurations in the JBoss server. The caNanoLab web application relies on a modified version of NCICBIIT CSM (Common Security Module) version 4.1 for user authentication and role-based authorization. For more information, see the CSM documentation at <a href="http://gforge.nci.nih.gov/frs/?group_id=12">http://gforge.nci.nih.gov/frs/?group_id=12</a>.</p> <p><b>Note:</b> You may skip this step if you have JBoss set up already from release 1.5.</p>

Step	Action
6	<p>Execute the Ant build script <code>build.xml</code> located at <code>&lt;CANANOLAB_SOURCE&gt;</code> with the target <code>reset-csm-db-connection</code>.</p> <p><i>Example:</i> Issue the following commands to execute the Ant script:</p> <pre>C:\&gt;cd C:\caNanoLab_1.5.1 C:\caNanoLab_1.5&gt;ant reset-csm-db-connection</pre> <p>Successful execution of the Ant script sets up the required CSM 4.1 configurations in the database.</p>
7	<p>When deploying both the caNanoLab application and the caNanoLab grid service (see page 12) to the same JBoss server, we recommend increasing the JBoss JVM heap size to 1G bytes and permanent generation (permgen) memory space to 512M bytes by updating the file <code>&lt;JBOSS_HOME&gt;/bin/run.conf</code> with the following <code>JAVA_OPTS</code>:</p> <pre>JAVA_OPTS="-Xms512m -Xmx1024m -XX:PermSize=256m - XX:MaxPermSize=512m -Dsun.rmi.dgc.client.gcInterval=3600000 - Dsun.rmi.dgc.server.gcInterval=3600000"</pre> <p>This configuration file is read when JBoss server starts.</p>
8	<p>When deploying the caNanolab application and the grid service in a production environment, we also recommend updating the default logging behavior of the JBoss server by turning off the unnecessary loggings thus reducing file system requirements for server log files.</p> <p>Please locate the file <code>&lt;JBOSS_HOME&gt;/server/default/conf/log4j.xml</code> and add to the beginning of the section <code>Limit categories</code> the following text:</p> <pre>&lt;!-- ===== --&gt; &lt;!-- Limit categories --&gt; &lt;!-- ===== --&gt;  &lt;!-- edited for caNanoLab --&gt; &lt;category name="org.hibernate"&gt;   &lt;priority value="WARN"/&gt; &lt;/category&gt; &lt;category name="org.springframework"&gt;   &lt;priority value="WARN"/&gt; &lt;/category&gt; &lt;category name="org.globus"&gt;   &lt;priority value="WARN"/&gt; &lt;/category&gt; &lt;category name="gov.nih.nci"&gt;   &lt;priority value="ERROR"/&gt; &lt;/category&gt; &lt;category name="org.displaytag"&gt;   &lt;priority value="ERROR"/&gt; &lt;/category&gt;</pre>


Step	Action
	<pre> &lt;category name="org.directwebremoting"&gt;   &lt;priority value="ERROR"/&gt; &lt;/category&gt; &lt;category name="net.sf.ehcache"&gt;   &lt;priority value="INFO"/&gt; &lt;/category&gt; &lt;category name="com.mchange.v2"&gt;   &lt;priority value="INFO"/&gt; &lt;/category&gt; &lt;category name="org.jboss"&gt;   &lt;priority value="INFO"/&gt; &lt;/category&gt; &lt;category name="CSM"&gt;   &lt;priority value="ERROR"/&gt; &lt;/category&gt; &lt;category name="Authorization"&gt;   &lt;priority value="ERROR"/&gt; &lt;/category&gt; &lt;!--end of edit for caNanoLab --&gt; </pre>

## Verification

Once the deployment war files have been deployed and the JBoss application server is correctly configured, you can now start the JBoss application server, which in turn starts the caNanoLab application.

Open the URL [http://<APP\\_SERVER\\_URL>/caNanoLab/](http://<APP_SERVER_URL>/caNanoLab/) (for example, <http://cananolab-dev.nci.nih.gov:19080/caNanoLab>). You should see a disclaimer page. Click on **CLICKING HERE** to go to the Welcome/Login page.


## User Provisioning Tool (UPT)

 <p><b>Creating User Accounts</b></p>	<p>Before users can log in to the caNanoLab application to submit and search data, you must first create their user accounts and assign them to the pre-defined user groups with the pre-defined roles to access the pre-defined protection groups. The caNanoLab application makes use of the NCICBIIT's User Provisioning Tool (UPT), a separate web application, for user account management. The concepts of users, groups, roles, protection groups are defined according the CSM/UPT principles. See the CSM documentation at <a href="http://gforge.nci.nih.gov/frs/?group_id=12">http://gforge.nci.nih.gov/frs/?group_id=12</a> for details on these concepts and the use of the UPT tool.</p>
--------------------------------------------------------------------------------------------------------------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

In release 1.5.1, as a part of the database seed data, a default user group `Public` has been created and has been assigned role `R` to protection group `protocol`, `sample` and `publication`, i.e. assigned access to search public protocols, samples and publications.

During the caNanoLab application start up, two more default users groups are auto-generated: `<application.owner>_Curator`, `<application.owner>_Researcher`, where the value for `<application.owner>` is specified in the file `build.properties` prior to building the application as described on page 6. At start up, `<application.owner>_Curator` is assigned role `CURD` to protection group `protocol`, `sample`, and `publication`, i.e. assigned access to submit and update protocols, samples and publications.

When a user in the `<application.owner>_Curator` group submits data (protocols, samples, publications) into the caNanoLab system, the data is automatically assigned to be visible for the two default groups.

<p><b>NOTE:</b></p> 	<p>In release 1.5.1, a user must be assigned as a caNanoLab administrator to see the <b>ADMINISTRATION</b> menu item in the application to log into the UPT tool or do update site preferences such as site logo.</p>
-------------------------------------------------------------------------------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

### UPT Example

The following steps illustrate an example use of the UPT tool to create a new user, to assign the user to be a caNanoLab administrator, and to assign the user to a pre-defined user group.


Step	Action
1	<p>Launch the UPT tool at <a href="http://&lt;APP_SERVER_URL&gt;/upt">http://&lt;APP_SERVER_URL&gt;/upt</a> and login as <b>superadmin/superadmin</b>. Use <b>csmupt</b> as the application name when prompted at the UPT log in.</p> <p><b>Note:</b> The user <b>superadmin</b> with initial password <b>superadmin</b> was created as a part of the database setup. Only <b>superadmin</b> can assign users to be caNanoLab administrators.</p>
2	<p>Once logged into the UPT tool, follow these steps to reset the password for <b>superadmin</b>:</p> <ol style="list-style-type: none"> <li>Select <b>User &gt; Select an Existing User</b>.</li> <li>Click on <b>Search</b> and select <b>superadmin</b> and click on <b>View Details</b>.</li> <li>Update the <b>User Password</b> field and <b>Confirm Password</b> field and click on <b>Update</b> to commit the change.</li> </ol>

Step	Action
3	<p>Logged in as the super admin, follow these steps to create a new user:</p> <ul style="list-style-type: none"> <li>d. Select <b>User &gt; Create a New User</b>.</li> <li>e. Create a new user account named <b>admin</b> with an initial password.</li> <li>f. Select <b>Application &gt; Select an Existing Application</b>; click <b>Search</b>.</li> <li>g. Select <b>caNanoLab</b> from the application list.</li> <li>h. Click <b>View Details</b>, then <b>Associated Admins</b>.</li> <li>i. Assign this user to be an administrator for the caNanoLab application.</li> <li>j. Click <b>Update Association</b> to commit the change.</li> </ul>
4	<p>Log out of the UPT tool and log back in as <b>admin</b>. Use <b>caNanoLab</b> as the application name when prompted at the UPT log in.</p>
5	<p>Select <b>User &gt; Select an Existing User</b>, and click <b>Search</b>. Select <b>admin</b> from the User list.</p> <ul style="list-style-type: none"> <li>a. Click <b>View Details</b>, then <b>Associated Groups</b>.</li> <li>b. Select <b>&lt;application.owner&gt;_Curator</b> from the pre-defined group list and assign it to the user.</li> <li>c. Click <b>Update Associations</b> to commit the change.</li> </ul>

---

Follow similar steps to create other application user accounts and to assign them to different users groups, as appropriate.

---

<p><b>NOTES:</b></p> 	<ul style="list-style-type: none"> <li>• In release 1.5.1, <b>superadmin</b> or caNanoLab administrators assign the initial password when creating a new user account.</li> <li>• The UPT tool doesn't allow users to manage their own passwords. Users can only update their passwords within the caNanoLab application by clicking on the <b>Update Password</b> link on the home page.</li> <li>• Only users assigned as caNanoLab administrators and assigned to the <code>&lt;application.owner&gt;_Curator</code> group can delete information from caNanoLab.</li> <li>• Publicly available data can be browsed through <b>Browse caNanoLab</b> links on the home page without logins.</li> <li>• Users logged into the application have to be assigned to one of the predefined user groups to be able access protected data.</li> <li>• A user can be assigned to multiple user groups. Roles and protection groups are assigned only to user groups, not users.</li> <li>• During sample creation, the organization name of the primary point of contact associated with the sample is automatically created as a new user group, and the new sample is automatically assigned role <b>R</b> to the newly created user group, i.e. any users later assigned to this new user group (through the UPT tool) would have read access to this sample.</li> </ul>
--------------------------------------------------------------------------------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

---

For more information about how to use the UPT tool for managing user accounts, contact NCICBIIT Application Support at [NCICBIIT@pop.nci.nih.gov](mailto:NCICBIIT@pop.nci.nih.gov) and request that the caNanoLab technical team give you a demonstration of the UPT tool in the context of the caNanoLab application.

---

## caNanoLab Grid Data Service

---

### Assumptions and Requirements

The JBoss server has been downloaded, installed and configured correctly, and the release 1.5.1 caNanoLab application has been correctly installed. If you had started the JBoss server to test the caNanoLab application installation, please shut it down now. In release 1.5.1, the caNanoLab grid service is configured to be deployed to the Tomcat embedded in the same JBoss server used for deploying the caNanoLab application. Users are no longer required to install a separate Tomcat server or to download and install Globus Toolkit.

In order for the grid service to successfully register with the NCICBIIT production index server, the server and the port number hosting the JBoss server must be publicly accessible.

---

### Installation and Deployment

Follow these steps to create and deploy a caNanoLab grid data service.

Step	Action
1	<p>If not done already, set up an environment variable <code>JBOSS_HOME</code> to point to the JBoss server install directory as in step 1 on page 6.</p> <p>If not done already, edit the Ant build properties file <code>build.properties</code> at <code>&lt;CANANOLAB_SOURCE&gt;</code> as in step 2 on page 6.</p> <p>If not done already, edit the Ant build properties file <code>deploy.properties</code> at <code>&lt;CANANOLAB_SOURCE&gt;</code> as in step 3 on page 7.</p> <p><b>Note:</b> the value of the property <code>application.owner</code> is used to construct the site name shown in the <b>Site</b> drop-down list on the caNanoLab home page.</p>

2	<p>Edit the Ant build properties file <code>grid.properties</code> at <code>&lt;CANANOLAB_SOURCE&gt;</code> by specifying values for the following properties:</p> <ol style="list-style-type: none"> <li><code>grid.indexserver</code>: the grid index server from which the caNanoLab grid data services can be automatically discovered, for example, <a href="http://cagrid-index.nci.nih.gov:8080/wsrf/services/DefaultIndexService">http://cagrid-index.nci.nih.gov:8080/wsrf/services/DefaultIndexService</a>, the NCICBIIT's production caGrid index server.</li> <li><code>grid.service.hostname</code>: virtual host name for the server hosting JBoss, for example, <code>cananolab.nci.gov</code></li> <li><code>grid.service.port</code>: port number configured for the Tomcat service embedded in JBoss, for example, <code>80</code></li> <li><code>country</code>: the country in which the caNanoLab grid service is running, for example, <code>US</code>.</li> <li><code>city</code>: the city at which the caNanoLab grid service is running, for example, <code>Rockville</code>.</li> <li><code>zipcode</code>: the zip code at which the caNanoLab grid service is running, for example, <code>20852</code>.</li> <li><code>state</code>: the state at which the caNanoLab grid service is running, for example, <code>MD</code>.</li> <li><code>street1</code>: the first part of the street at which the caNanoLab grid service is running, for example, <code>2115 East Jefferson Street</code>.</li> <li><code>street2</code>: the second part of the street at which the caNanoLab grid service is running.</li> <li><code>affiliation</code>: the affiliation of the person who is responsible for maintaining the caNanoLab grid service, for example, <code>NCICBIIT</code>.</li> <li><code>email</code>: the email of the person who is responsible for maintaining the caNanoLab grid service, for example, <a href="mailto:NCICBIIT@pop.nih.gov">NCICBIIT@pop.nih.gov</a></li> <li><code>phone</code>: the phone number of the person who is responsible for maintaining the caNanoLab grid service.</li> <li><code>first.name</code>: the first name of the person who is responsible for maintaining the caNanoLab grid service.</li> <li><code>last.name</code>: the last name of the person who is responsible for maintaining the caNanoLab grid service.</li> </ol> <p><b>Note:</b> Properties c through g are grid service metadata that are required for the grid index server to visually map the location of a grid service in the caGrid portal, and are useful during grid service auto-discovery. Properties h through m are optional, but useful information for others to look up the grid service.</p>
---	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------



3	<p>Execute the Ant build file <code>build.xml</code> located at <code>&lt;CANANOLAB_SOURCE&gt;</code> with the target <code>deploy-grid</code>. For example,</p> <pre>C:\&gt;cd C:\caNanoLab_1.5.1 C:\caNanoLab_1.5.1&gt;ant deploy-grid</pre> <p>Successful execution of the Ant build script generates a war file <code>wsrf-canano.war</code> at the <code>&lt;CANANOLAB_GRID_SOURCE&gt;/build</code> directory and places the war file to the JBoss deploy directory, for example, <code>&lt;JBOSS_HOME&gt;/server/default/deploy</code>.</p>
4	<p>Execute the Ant build script <code>build.xml</code> located at <code>&lt;CANANOLAB_SOURCE&gt;</code> with the target <code>prepare-globus-lib</code>. For example,</p> <pre>C:\&gt;cd C:\caNanoLab_1.5.1 C:\caNanoLab_1.5.1&gt;ant prepare-globus-lib</pre> <p>Successful execution of the Ant build script extracts eight jar files for deploying Globus to the JBoss server at <code>&lt;JBOSS_HOME&gt;/server/default/lib</code>.</p>
5	<p>If not done already, follow step 7 and 8 on pages 9 and 10 to configure appropriate memory and logging settings for the JBoss server.</p>

## Verification

Once the grid service war file has been deployed and you have correctly configured JBoss server with required Globus libraries, you can now start the JBoss server, which in turn starts the caNanoLab application and the caNanoLab grid service. To verify the creation of the caNanoLab grid data service, open the URL [http://<APP\\_SERVER\\_URL>/wsrf-canano/services/cagrid/CaNanoLabService](http://<APP_SERVER_URL>/wsrf-canano/services/cagrid/CaNanoLabService), for example, <http://cananolab-dev.nci.nih.gov:19080/wsrf-canano/services/cagrid/CaNanoLabService>. You should see a page that says:

```
cagrid/CaNanoLabService
```

```
Hi there, this is an AXIS service!
```

It takes more than ten minutes for the grid service to successfully register in the index server. Wait some time and go into the caGrid portal at <http://cagrid-portal.nci.nih.gov/web/guest/home>. You should see your grid service showing up on the portal map. This completes the release 1.5.1 installation.

**Note:** After successful registration with the production grid index server, your grid service name will show up as a remote site in the **Site** drop-down list on the home pages of other caNanoLab deployments.

For more information and help on caGrid service configuration and deployment, refer to the caGrid wiki page: <http://www.cagrid.org>.

---

If you need further assistance on setting a caNanoLab grid service, contact NCICBIIT Application Support at [NCICBIIT@pop.nci.nih.gov](mailto:NCICBIIT@pop.nci.nih.gov) and request that the caNanoLab technical team help you with caNanoLab grid service.

---

## Contacting Application Support

---

**NCICBIIT  
Application  
Support**

<http://NCICBIIT.nci.nih.gov/NCICBIIT/support>  
[NCICBIIT@pop.nci.nih.gov](mailto:NCICBIIT@pop.nci.nih.gov)  
Telephone: 301-451-4384  
Toll free: 888-478-4423

---